



Fast scanning techniques for cancer therapy with hadrons

—

a domain of cyclotrons

Marco Schippers,
David Meer, Eros Pedroni

(thanks to all colleagues in the proton therapy field)





PAUL SCHERRER INSTITUT

accelerator

energy selection

beam transport

gantry / fixed hor. line




To understand the need for a cyclotron.....

- Dose delivery techniques
- Fast 3D Pencil beam scanning
- Pencil beam scanning results at PSI
- Consequences for accelerator

Dose delivery techniques



Dose delivery techniques: Depth



Tumor **distal edge**

→ Range

→ Maximum Energy
per field → „slow“ (sec)

Tumor **thickness**

→ spread-out Bragg peak

→ energy modulation

During trmt → „fast“ (<0.1 sec)

Methods:

- 1) at accelerator
- 2) just before patient (in “nozzle”)

Dose delivery techniques: Depth

Vary energy at accelerator:


FFAG, Linear accelerators:

- Pulsed machines
- Energy adjustable **per pulse**



Suitable for E-modulation

if Energy and Dose are correct at each pulse

Dielectric Wall Acc




Cyclotron+Linac BOoster




Dose delivery techniques: Depth

Vary energy at accelerator: **synchrotron**

Vary energy **at** each spill:



Vary energy **during** a spill:






Iwata et al., IPAC'10, Kyoto, Japan MOPEA008

→ Sets range

→ **Additionally**: Fast energy modulation in nozzle

Vary energy at accelerator: cyclotron



Rolled-up wedge degrader
220-70 MeV (IBA)

- Sets range
- **Additionally:** Fast energy modulation in nozzle


Dose delivery techniques: Depth

Vary energy at accelerator: cyclotron

PSI and VARIAN: fast degrader + laminated magnets



Multiple wedge degrader
238-70 MeV





→ 5 mm Δ Range in 50 ms

Gantry-2: No Additional energy modulation in nozzle


Dose delivery techniques: Depth

PSI fast degrader system: in depth 5 mm in 80 ms

- fast treatment
- fast treatment area switching



Fast magnets: quick area switching




Dose delivery techniques: Depth


Degrader purpose: **decrease energy**

however: - **energy spread increases**


degrader system



Collimators define
transmitted beam size



- **beam size increases** due to multiple scattering
- **beam loss** due to nuclear reactions in degrader




Van Goethem et al., Phys. Med. Biol. 54 (2009)5831


Dose delivery techniques: Depth

Degrader purpose: **decrease** energy

however: - **energy spread** increases




degrader system



Collimators define
transmitted beam size


- **beam size** increases due to multiple scattering
- **beam loss** due to nuclear reactions in degrader

→ Beam intensity from cyclotron must be high enough



Van Goethem et al., Phys. Med. Biol. 54 (2009)5831

Dose delivery techniques: Depth



Degrader and collimator: careful material choice



=> Activation decays quickly

=> Amount of activated material is limited

Dose delivery techniques: Depth

Vary energy in nozzle (cyclotron and synchrotron)

Max penetration depth: bolus or compensator





But: material in front of patient
- increases scatter → unsharp edges
- enter treatment room to change parts

Dose delivery techniques: Depth

Vary energy in nozzle (cyclotron and synchrotron)


Energy modulation: rotating wheel or insertable plates




But: material in front of patient
 - increases scatter → unsharp edges
 (- enter treatment room to change parts)

Dose delivery techniques: Depth

Energy modulation **in nozzle**: no beam analysis




Energy modulation **upstream**: includes beam analysis



Dose delivery techniques: lateral


Scatter technique




Pencil beam scanning



Scattered beam




Spot scanning



Dose distribution of
scattered beam

Spot scanning: step & shoot



Beam size 7 mm FWHM
5 mm steps
10'000 spots/liter (21 x 21 x 21)
Dose painted only once
~1 Gy / liter/minute



Requirements for accelerator:

- stable beam position


Scattering – Scanning

Scattered beam


Scanned beam



Intensity Modulated Proton Therapy




Intensity Modulated Proton Therapy




calculation made
by T.Lomax, PSI

Intensity Modulated Proton Therapy




calculation made
by T.Lomax, PSI

Intensity Modulated Proton Therapy




calculation made
by T.Lomax, PSI

Intensity Modulated Proton Therapy



calculation made
by T.Lomax, PSI

Intensity Modulated Proton Therapy




Total dose
of all 4
beam directions

calculation made
by T.Lomax, PSI

The problem in dynamical treatments:


Organ movement

➔ Danger to underdose and overdose



The problem in dynamical treatments:

Organ movement




→ Danger to underdose and overdose

Solutions:

The problem in dynamical treatments:

Organ movement



➔ Danger to underdose and overdose


Solutions:

- Beam gating
- Multiple scans of tumour
- „follow“ tumour

⇒ increase scan speed laterally + in depth

The problem in dynamical treatments:

Organ movement



➔ Danger to underdose and overdose

Solutions:



- Beam gating
- Multiple scans of tumour
- „follow“ tumour

⇒ increase scan speed laterally + in depth

Fast pencil beam scanning in 3D

Cont. scanning "TV" mode





kHz-Intensity modulation



- **7 s** for a **1 liter** volume.
- Target **repainting**:
- 15-30 scans / 2 min.


Intensity control

Max. intensity set by:
proton source






Deflector plate:
sets requested intensity

- within 50 μ s
- 5% accuracy





+ slits

Intensity (Vdefl) → intensity modulation



Intensity correction for Transmission(E)

Fast intensity regulation over factor 20 needed.




"quick and dirty"
brute force

Project coming up: Fast change of I_{max}


Regulation of beam intensity over **factor 20** needed.

Idea: **vary HF power** in stead of Qpole defocusing




Study started to increase window of high extraction efficiency


PSI Gantry-2: fast 3D scanning



$\sigma = 3-4.5$ mm over 230-100 MeV




Fast rescanning demo:
500 planes of 6x8 cm in depth < 1 min.



scanning magnets (x & y)
for 2D // scanning

R=3.5 m



PSI "Gantry-2"

Eros Pedroni

David Meer

Requirements for accelerator

Only one conclusion....



Fast pencil beam scanning

- reliable, robust, easy to tune
- enough intensity
- stable beam position
- CW beam
- stable beam intensity
- fast control of beam intensity
- fast control of beam energy
- small pencil beam
- controlled by user (tune library)

Cyclotrons

You can't get around it