



First Measurements of Photoelectron Transverse Energy Distribution Curve using TESS

D.V. Gorshkov, H.E. Scheibler, A.S. Terekhov

Rzhanov Institute of Semiconductor Physics SB RAS, Novosibirsk, Russia

R.J. Cash, B.D. Fell, L.B. Jones, K.J. Middleman,

B.L. Militsyn, T.C.Q. Noakes

ASTeC, STFC Daresbury Laboratory, Daresbury, United Kingdom

TESS - Transverse Energy Spread Spectrometer, located at Daresbury Laboratory, UK.









ERL-2013











Electron transit time from photocathode to MCP front surface τ is equal to

$$\tau = d \times \sqrt{\frac{2 m_{e}}{e U_{acc}}} \times \left(\sqrt{1 + \frac{\varepsilon_{lon}}{e U_{acc}}} - \sqrt{\frac{\varepsilon_{lon}}{e U_{acc}}}\right)$$

where d – distance between photocathode and MCP, U_{acc} – accelerating voltage, applied between photocathode and MCP, \mathcal{E}_{lon} – longitudinal energy of the electron, m_e- free electron mass, e – electron charge.

In photocathode plane electron moves a distance r

$$r = \sqrt{\frac{\varepsilon_{tr}}{2 m_e}} \times \tau$$
 and $\varepsilon_{tr} = \frac{m_e r^2}{2 \tau^2}$





Conversion of $I_{el}(r)$ to $N(\epsilon_{tr})$

Amount of electrons which hit the ring with radius r and thickness ∂r is equal to

 $\partial N = A \times 2 \pi r \, \partial r \times I_{el}(r)$

$$N(\varepsilon_{tr}) = \frac{\partial N}{\partial \varepsilon_{tr}} = \frac{A 2 \pi r \partial r I_{el}(r)}{\frac{m_e r}{\tau^2} \times \partial r} \sim I_{el} \left(r = \sqrt{\frac{\varepsilon_{tr} 2 \tau^2}{m_e}} \right)$$

If we can neglect ε_{lon} ($\varepsilon_{lon} << e U_{acc}$)

$$\tau = d \times \sqrt{\frac{2 m_e}{e U_{acc}}}$$
 and $N(\varepsilon_{tr}) \sim I_{el} \left(r = 2 d \times \sqrt{\frac{\varepsilon_{tr}}{e U_{acc}}}\right)$

TESS – The Transverse Energy Spread Spectrometer



- Independent control of cathode & grid potentials
- Ability to cool the photocathode ٠ to LN₂ temperature
- Ability to change the ٠ wavelength of the laser beam
- Ability to 'poison' the cathode ٠ through a leak valve

Phosphor screen

Retarding grids











TEDC measurements at different accelerating voltages



17.06.2013 p-GaAs(Cs,O) photocathode



ERL-2013



TEDC measurements at different accelerating voltages



17.06.2013 p-GaAs(Cs,O) photocathode

Screen luminescence profiles





TEDC measurements at different accelerating voltages



17.06.2013 p-GaAs(Cs,O) photocathode

Transverse energy distribution curves



 $d=33\ mm\ \lambda=635\ nm\ T=300\ K$





TEDC measurements at different drift distances



20.06.2013 p-GaAs(Cs,O) photocathode



$$U_{acc} = 205 \text{ V} \ \lambda = 635 \text{ nm} \ T = 300 \text{ K}$$





TEDC measurements at different drift distances



20.06.2013 p-GaAs(Cs,O) photocathode

Screen luminescence profiles





TEDC measurements at different drift distances



20.06.2013 p-GaAs(Cs,O) photocathode

Transverse energy distribution curves



 $U_{acc} = 205 \text{ V} \quad \lambda = 635 \text{ nm} \quad T = 300 \text{ K}$





TEDC measurements with different light wavelength



21.06.2013 p-GaAs(Cs,O) photocathode $d \approx 43 \text{ mm } T = 300 \text{ K}$ $\lambda = 532 \text{ nm}$ $\lambda = 635 \text{ nm}$



$$U_{acc} = 60 V$$





TEDC measurements with different light wavelength



21.06.2013 p-GaAs(Cs,O) photocathode $d \approx 43 \text{ mm } T = 300 \text{ K}$



MTE₆₃₅= 45 ± 7 meV MTE₅₃₂= 100 ± 15 meV

ERL-2013



Summary



- Transverse Energy Spread Spectrometer (TESS) A new set-up for measuring transverse energy distribution curves (TEDC) of electrons emitted by a reflectance-mode photocathode, has been developed, constructed and commissioned.
- TESS enables one to measure:
 - TEDC of the electrons emitted by all kinds of solid-state reflectancemode photocathodes: semiconductor photocathodes, metal photocathodes, bi-alkali and multi-alkali photocathodes;
 - TEDC dependence on the wavelength of exciting radiation;
 - TEDC during photocathode degradation, including poisoning of the photocathode by the gases, precisely introduced into the vacuum chamber;
 - TEDC of the electrons, emitted from photocathodes at RT and LN₂ temperatures.

