

37th International Free Electron Laser Conference

Direct observation of bond formation in solution with femtosecond X-ray scattering

2015. 8. 28

Kim, Kyung Hwan

Center for Nanomaterials and Chemical Reactions,
Institute for Basic Science (IBS)

Contents


1. Introduction

- TRXSS method – How to measure structural changes in solution

2. Application

- Tight Au-Au bond formation in solution


Scattering from a molecule





- **X-ray scattering** is one of the most useful technique for **revealing structure** of sample.

Scattering from solution


- Too many atomic pairs.
- Impossible to analyze.




$[\text{Au}(\text{CN})_2^-]_3$ in water




Pump-probe scheme




Pump-probe scheme




Pump-probe scheme



Pump-probe scheme





Pump-probe scheme




Signals only from the structural changes
All (known or unknown) backgrounds are subtracted.

Example: I₂

The **bondlength** of I₂ molecule (2.67 Å) is measured with **high accuracy**.




F. T.
→




Example: I₂

$$q\Delta S(q)_{theory} = 2S(q)_I - S(q)_{I_2} = F_I(q)^2 \frac{\sin(qr)}{qr} \quad (r = 2.67 \text{ \AA})$$




$$r\Delta R(r, t) = \frac{1}{2\pi^2} \int_0^\infty q\Delta S(q, t) \sin(qr) \exp(-q^2\alpha) dq$$



F. T.
→



Multiple number of atom pairs



Synchrotron Sources

ESRF



APS



KEK




- Location : Grenoble, France
- Beamline : ID09B
- X-ray : 100ps, 18keV, 1.1×10^9 photon/pulse
- Coworker : M. Wulff


- Location : Argon, USA
- Beamline : 14-ID
- X-ray : 100ps, 12keV, $\sim 10^9$ photon/pulse
- Coworker : K. Moffat

- Location : Tskuba, Japan
- Beamline : NW14
- X-ray : 100ps, 15keV, $\sim 10^9$ photon/pulse
- Coworker : S. Adachi


Applications



JACS, 130, 5834-5835 (2008)




Science, 309, 1223-1227 (2005)
JPCA 1209 10451-10458 (2005)
JCP, 124, 124504 (2006)




ACIE, 47, 5550-5553 (2008)
JACS, 132, 2600-2607 (2010)



JACS, 129, 13584-13591 (2007)




JCP, 124, 034501 (2006)
PCCP, 17, 8633-8637 (2015)



ACIE, 47, 1047-1050 (2008)



PNAS, 103, 9410-9415 (2006)



PRL, 110, 165505 (2013)
ChemPhysChem, 14, 3687-3697 (2013)

Femtosecond TRXSS using X-ray free electron laser


	3 rd generation	XFEL
Temporal duration	100 ps	~ 0.1 ps
Photons/pulse	< 10 ⁹	~ 10 ¹²

SACLA, JAPAN



- Real-time tracking of bond-making process.

Timescale of ultrafast sciences




Elementary steps of chemical reactions

Unimolecular & Bimolecular reaction



Real-time tracking of "bond-breaking" processes in solution





Ihee *et al.*, *Science*, **309**, 1223 (2005).

- Bond-breaking processes **can be synchronized** with laser excitation.
- Relatively **easy to follow** with TR techniques.


Real-time tracking of "bond-making" processes in solution

$h\nu$




- Bond-making process >> Diffusion rate
- Hard to monitor ultrafast bond-making process with TR techniques.

Real-time tracking of "bond-making" processes in solution




- Bond-making processes can be synchronized with laser excitation if the reaction parties are prepared in the same cage. ($[\text{Au}(\text{CN})_2^-]^3$ system)

Auophilic interaction of Gold




- Special van der Waals interaction induced by relativistic effect.
- Longer Au-Au distance (3.0 ~ 3.6 Å) than covalent bond (~2.7 Å).



Ganglong Cui *et al.*, *Angew. Chem. Int. Ed.*, **52**, 10281 –10285 (2013).



Laser induced bond formation



Ganglong Cui *et. al.*, *Angew. Chem. Int. Ed.*, **52**, 10281 –10285 (2013).

Lack of structural sensitivity of previous study

Controversy over transient structure




- 0.5 ps, 2.1 ps, 2ns kinetics
- Bent-to-linear relaxation with 2.1 ps.

Munetaka Iwamura *et al.*, *J. Am. Chem. Soc.* **135**, 538 (2013)

- Bent-to-linear relaxation occur within 500 fs.

Ganglong Cui *et al.*, *Angew. Chem. Int. Ed.*, **52**, 10281 (2013).

Experimental scheme



- Typical TRXL experimental scheme was used.
- We covered from -800 fs to $1 \mu\text{s}$ time range. (SACLA + KEK)

Experimental Conditions

Sample

- ◆ 300 mM $\text{Au}(\text{CN})_2^-$ in water.


X-ray

- ◆ 15 keV
- ◆ Pulse duration: <100 fs (SACLA), ~ 100 ps (KEK).


Laser

- ◆ 267 nm
- ◆ Pulse duration: ~100 fs

XFEL + Synchrotron scheme




TRXL signal from -800 fs to 1 μ s




- Clear oscillatory features which varies with time.
- High S/N.



TRXL data analysis scheme




Fourier Transformation
Theoretical calculation
1) solute-only term
2) solute-solvent cross term
3) solvent heating



Structure reconstruction




TRXL data analysis scheme





II




+



+




Time-dependent RDFs





- The peaks in RDFs directly represent the Au-Au distances.
- Extract kinetics and species-associated RDFs.

Concentration changes





- 1.6 ps (S_1 to T_1), 3 ns (T_1 to tetramer), 100 ns (tetramer to S_0)
- Timescales are well matched with previous studies.

Species-associated RDFs: S_0





- S_0 :
Loosely bounded (~3.6 Å), bent structure ($R_{12}+R_{23} > R_{13}$).

Species-associated RDFs: S_1





- S_1 : covalent bond is formed ($\sim 2.8 \text{ \AA}$), linear structure ($R_{12}+R_{23}=R_{13}$), formed within $\sim 500 \text{ fs}$

Species-associated RDFs: T₁



- T₁:
shorter bond lengths
linear structure
($R_{12}+R_{23}=R_{13}$),
formed with a time constant of **1.6 ps**



Species-associated RDFs: tetramer



- Tetramer:
addition of a monomer,
formed with a time
constant of **3 ns**

Lack of structural sensitivity of previous study

Controversy over transient structure




- 0.5 ps, 2.1 ps, 2ns kinetics
- Bent-to-linear relaxation with 2.1 ps.

Munetaka Iwamura *et al.*, *J. Am. Chem. Soc.* **135**, 538 (2013)

- Bent-to-linear relaxation occur within 500 fs.

Ganglong Cui *et al.*, *Angew. Chem. Int. Ed.*, **52**, 10281 (2013).

Bent? Linear?




- Linear structure fits the experimental curve much better both in r and q space.

Summary: Molecular movie



K.H. Kim *et al.*, *Nature*, **518**, 385 (2015).

Summary



Structural parameters

Species	R_{12}	R_{23}	R_{13}	R_{34}
S_0	3.87 (± 0.04) Å	3.30 (± 0.06) Å	5.56 (± 0.11) Å	--
S_1	2.82 (± 0.04) Å	2.81 (± 0.03) Å	5.63 (± 0.09) Å	--
T_1	2.71 (± 0.03) Å	2.70 (± 0.05) Å	5.41 (± 0.11) Å	--
Tetramer	2.89 (± 0.06) Å	2.62 (± 0.06) Å	--	2.88 (± 0.04) Å

Conclusion

- We demonstrated the capability of femtosecond TRXL by elucidating the overall mechanism for the formation of Au-Au covalent bonds in the $[\text{Au}(\text{CN})_2^-]^3$ complex with rich structural information.
- Femtosecond TRXL offers an opportunity of visualizing the entire process of photoinduced reactions in real time and real space.

Acknowledgement

Center for Nanomaterials and Chemical Reactions

Prof. Hyotcherl Ihee
Jong Goo Kim
Key Young Oang
Tae Wu Kim
Hosung Ki
Junbeom Jo
Sungjun Park

KEK

Prof. Shin-ichi Adachi
Prof. Shunsuke Nozawa
Dr. Tokushi Sato

Inha University
Prof. Jeongho Kim
The Catholic University
Prof. Joonghan Kim

Funding

This work was supported by IBS-R004-G2

SACLA

Prof. Tetsuya Ishikawa
Prof. Changyong Song
Dr. Makina Yabashi
Dr. Tokushi Sato
Dr. Takahiro Sato
Dr. Kanade Ogawa
Dr. Tadashi Togashi
Dr. Kensuke Tono



Thank you.