

Extremely Low Emittance Beam Size Diagnostics with Sub-Micrometer Resolution Using Optical Transition Radiation

K. Kruchinin, P. Karataev, S.T. Boogert, L. Nevay

John Adams Institute at Royal Holloway, Egham, Surrey, TW20 0EX

A. Aryshev, M. Shevelev, N. Terunuma, J. Urakawa

KEK, 1-1 Oho, Tsukuba, Ibaraki 305-0801, Japan

B. Bolzon, T. Lefevre, S. Mazzoni


CERN, Geneva, Switzerland



Outlines


- Introduction
- Experimental setup
- Data analysis and calibration
- Monitor tuning
- Future improvements
- Summary

Introduction





- Transition radiation (TR) appears when a charged particle crosses a boundary between two media with different dielectric constants.
- The resolution is determined by the source dimensions induced by a single particle plus distortion caused by the optical system (diffraction of OTR tails)


ATF II Extraction Line



Experimental Setup



OTR screen



- Lens - "CVI Laser Optics" cemented achromat, $f=120\text{mm}$, $\phi=30\text{mm}$
- CCD Camera - SBIG-ST8300M with $5.4 \mu\text{m}$ pixel size, 3352×2532 pixel array and $\sim 50\%$ quantum efficiency

OTR Image






Image of the OTR spot taken with linear polarizer and 550 ± 20 nm optical filter

Horizontal Projection



Vertical Projection



PSF-like fit function:

$$f(x) = a_0 + \frac{a_1 \left(a_4 + (x - a_3)^2 \right)}{1 + (a_2 (x - a_3))^4}$$

Fit parameters:

- a_0 is the vertical offset of the distribution with respect to zero
- a_1 is the amplitude of the distribution
- a_2 is the smoothing parameter
- a_3 is the horizontal offset of the distribution with respect to zero
- a_4 is the distribution width

PSF-like Fit Function

Analytical calculation of:

- Minimum to maximum ratio

$$I_{min}/I_{max} = \frac{2a_2^2 a_4}{2a_2^2 a_4 + \sqrt{1 + a_2^4 a_2^2}}$$

- Distance between peaks

$$\frac{2\sqrt{-a_2^2 a_4 + \sqrt{1 + a_2^4 a_2^2}}}{a_2}$$

- Simpler expressions for error calculation

Self Calibration Procedure


- In the whole data set find a file with smallest I_{min}/I_{max}
- Regenerate fit curve $f(x)$ with errors for the calibration file substituting zeros for horizontal offset a_0 and smoothing parameter a_4
- Convolute fit with Gaussian distribution as follows:

$$F(x_j) = \frac{\sum_{i=1}^N f(x_i) \exp\left(-\frac{(x_j - x_i)^2}{2\sigma^2}\right)}{\sum_{i=1}^N \exp\left(-\frac{(x_j - x_i)^2}{2\sigma^2}\right)}$$


- Propagate errors through convolution.
- Repeat convolution N times varying σ from 0 to σ_m with a fine step.
- For each iteration. find I_{min}/I_{max} and calculate its errors resulting in calibration curve

Self Calibration Procedure

Beam size effect



Calibration curve




Fit function:


$$f(x) = a_0 + \frac{1}{a_1} \left(-\ln \left(1 - \frac{x}{a_2} \right) \right)^{a_3} + a_4 x^{12}$$

Best Quadrupole Scan

Minimum measured beam size was: $0.754 \pm 0.034 \mu\text{m}$



LW and OTR Comparison



More about LW see at [MOPF16](#)

Monitor Tuning and Optimization

Off - line:


- Alignment of the optical system and magnification factor measurements

On - line:


- Find OTR spot
- Optimization
 - Polarizer angle
 - Focusing scan
 - Image rotation scan
- Quadrupole scan
 - Find a file with the best visibility and perform calibration

Calibration of The Optical System


Image of the OTR screen



Calibration curve



Vertical projection




Fit function:

$$f(x) = a_0 + \frac{a_1}{1 + \exp\left(\frac{x-a_3}{a_2}\right)}$$


Measured magnification factor of
the system: 7.17

Focusing



11.70 mm 11.32 mm 11.54 mm 11.42 mm 11.30 mm



Focus-scan:



Optical Effects Study

Diffraction

Chromatic aberrations



Lens

"Sigma Koki" achromatic doublets:

DLB-30-120-PM

Properties

$f=120 \text{ mm}$, $\phi=30 \text{ mm}$

DLB-30-100-PM

$f=100 \text{ mm}$, $\phi=30 \text{ mm}$

"CVI Laser Optics" achromat lens

CVI LAO-120-30

$f=120 \text{ mm}$, $\phi=30 \text{ mm}$

"Sigma Koki" plano-convex lens:

SLB-30-100-PM



$f=100 \text{ mm}$, $\phi=30 \text{ mm}$

Future Improvements

- Using simulation tools such as ZEMAX in order to better understand the PSF behavior and optimize the optical system (see [MOPF04](#))
- Apply multi-element or reflective optics in order to reduce the resolution even further
- Upgrade experimental hardware (CCD camera, DAQ)
- Efforts towards automation (shot by shot beam size measurements)

ZEMAX Simulations (MOPF04)

- Zemax "Physical Optics Propagation" simulations: Propagation of the OTR source using diffraction laws (near field conditions) through real commercial lenses
- PSF simulation: Propagation of vertically polarized electric field for a single electron



Summary

- Collected Data shows a good progress in optimization of PSF-like OTR monitor system
- In order to improve the beam size measurement technique additional efforts toward the optimization of the optical system, and better understanding of the beam size effect has been taken.
- Work on analysis and simulations shows good agreement with experimental data