

The application

This system is designed for a synchrotron experimental station (a "beamline"). In this case, its purpose is SAXS/WAXS (Small- and Wide-angle X-ray scattering) measurements.


The system contains 2 main detectors: a Pilatus 300K (Dectris) to measure SAXS and a Pilatus 100K to measure WAXS. Pilatus 300K is configured to collect images at a given rate and sends a trigger TTL signal while an acquisition is being done.

A module of distributed I/O over EtherCAT, cRIO (NI), detects the trigger, opens the fast shutter and sends a TTL trigger to Pilatus 100K as it can start acquisition.

Two photo-diodes provide the beam intensity data through a SR 570 preamplifier that converts the current into voltage. cRIO reads this voltage and calculates its average while the trigger signal is high. The voltage average data is sent to PXI to be distributed to EPICS clients over the network.

Hardware used by LNLS SAXS1 beamline



cRIO with a FPGA algorithm: receive trigger signal → open shutter → collect voltage from SR 570 preamplifier until trigger signal is down → close shutter → calculate voltage average during the time the trigger signal was high → save the data



FPGA code running in cRIO


Sum of the voltage from Stanford until Pilatus 300K sets the trigger signal low. After this, the average is calculated.

Receive Scan Engine raising edge and transfers voltage average from FIFO to 100 shared-variables to EtherCAT.




Software in PXI

Real-time loop synchronized with Scan Engine receives data from EtherCAT and saves in an internal FIFO. This FIFO is read inside PXI and passed to a Linux IOC by a shared-memory.



Results

Trigger signal that is sent to the system, with 1.25 μ s high and 20 μ s low, resulting in a total period of 21.25 μ s.



Signal that indicates when FPGA stopped to read the voltage signal and detects that the trigger signal was low.

The minimum period we get in the worst case was **21.25 μ s**.