

IH-DTL as a compact injector for a heavy-ion medical synchrotron

Y. Iwata, S. Yamada, T. Murakami, H. Ogawa,
T. Fujisawa, N. Miyahara, K. Yamamoto,
S. Hojo, Y. Sakamoto, and M. Muramatsu
National Institute of Radiological Sciences (NIRS)

T. Fujimoto and T. Takeuchi
Accelerator Engineering Corporation (AEC)


T. Mitsumoto, H. Tsutsui, T. Ueda, and T. Watanabe
Sumitomo Heavy Industries Ltd. (SHI)

Contents

- Objective
- Design
- Beam acceleration tests
- Summary & conclusion

Objective

- Heavy Ion Medical Accelerator in Chiba (HIMAC)
- $E/A=800$ MeV for $q/m=1/2$
- Cancer treatments using energetic carbon ions
- 4,000+ patients (June, 1994 ~ present)



Widespread use of heavy-ion therapy

Present therapy complex

- Total size of complex: large
120 m X 60 m
- Construction costs: expensive
>\$300M

Widespread use of heavy-ion therapy

Developments of


■ Compact complex

55 m x 45 m

■ Compact injector

$L_{\text{tot}} < \sim 10 \text{ m}$

(Present: $> 32 \text{ m}$)



Required specifications

- Ion kind: Carbon ($^{12}\text{C}^{4+}$)
- Intensity: $>200 \text{ e}\mu\text{A}$
- Energy: $\sim 4 \text{ MeV/u}$
 - Injection energy of a synchrotron
 - Charge changing cross section
($^{12}\text{C}^{4+}$ to $^{12}\text{C}^{6+}$)

Design of compact injector

■ Ion source:

ECR Ion-Source (ECRIS)

- Technology is well-established
- Stable and reliable


Design of compact injector

■ Linacs:

RFQ linac

- Matching with the ECRIS
- Technology is well-established

Drift-Tube-Linac (DTL): APF IH-DTL




Compact ECRIS

■ Permanent magnets

- no P.S. for magnets
- no cooling water
- Easy maintenance
- Compact


■ Intensity of $^{12}\text{C}^{4+}$

>400 e μ A ($V_{\text{ext}} \sim 30$ kV)



RFQ linac


- 4-vane structure
- 200 MHz
- Beam energies
 - Injection: 10 keV/u
 - Extraction: 608 keV/u
- Cavity size
 - Length: 2.5 m
 - Diameter: 41 cm



APF IH-DTL

- **Beam focusing**
 - Alternating-Phase-Focusing (APF)
- **Cavity structure**
 - Interdigital H-mode (IH) structure
- **200 MHz**
- **Beam energies**
 - Injection: **608 keV/u**
 - Extraction: **4.0 MeV/u**

Principle of the APF method



Both stability of motion is obtained just with the rf acceleration field!

APF method

- Proposed in 50s
- Attractive features:
 - No focusing element
 - Simple cavity
 - Easy operation
- No practical APF linac exists
 - Weak focusing force
 - Design of synchronous phase

Phase array

■ Beam dynamics

- Array of synchronous phase

■ Simple case


- $\phi_s = +30, -30, +30, -30, \dots$ (FDFD)
- $\phi_s = +30, +30, -30, -30, \dots$ (FFDD)

■ Small acceptance

- $\beta_{\text{ext}}/\beta_{\text{inj}} \sim 2.6$
- Stable region changes with β


Phase array

- Sinusoidal-like function
- Iterative simulations
 - Tracking simulations
 - Gap voltages (3D field-calculations)
- Sufficient acceptance
 - Transmission~99.6%




IH structure

- Proposed in 50s
- Attractive features
 - High shunt impedance
(a few $\times 100 \text{ M}\Omega/\text{m}$)
 - Compact
- Electromagnetic field
 - 3D field solver



IH structure

- Weak focusing force of APF
 - $|\Delta V_g/V_g| < 2\%$
 - Precision of 3D solver
- Design of tuners
- Construction of model cavity





Model cavity

■ Full-scale model





Model cavity

- Full-scale model
- 3D field calculation
 - Micro Wave Studio (MWS)





Model cavity


- Full-scale model
- 3D field calculation
 - Micro Wave Studio (MWS)
- Field measurements
 - Perturbation method



Cavity of APF IH-DTL




Cavity of APF IH-DTL





Tuning of electric field

- **Field measurements**
 - Perturbation method
- **Designed field**
 - Calculated with MWS





Tuning of electric field

- **Field measurements**
 - Perturbation method
- **Designed field**
 - Calculated with MWS
- **Fine tuning**
 - **Good agreement**



High-power tests


- Required power was successfully fed into the cavity.



Features of APF IH-DTL

<Use of the APF method>

- **No focusing element**
 - Simple & cost-effective cavity
- **Small drift tubes**
 - Higher operating-frequency
 - Lower injection-energy
- **Easy operation**




<Use of IH cavity>

- **Compact cavity**
 - Length: **3.5 m**
 - Diameter: **0.44 m**
- **High shunt impedance**
 - Z_{eff} : **110 M Ω /m**
 - Required power: **~360kW**




Beam acceleration test (1)

■ ECR and RFQ linac




Beam acceleration test (1)

■ Transmission: 81%






Beam acceleration test (1)

■ Transmission: 81%



Beam acceleration test (1)



- Transmission: 81%
- Emittance: $1 \pi \cdot \text{mm} \cdot \text{mrad}$
(normalized 90% value)



Beam acceleration test (1)





■ Transmission: 81%

■ Emittance: 1.5 mm mrad




Beam acceleration test (1)

- Transmission: 81%
- Emittance: $1 \pi \cdot \text{mm} \cdot \text{mrad}$
(normalized 90% value)
- Energy: $\langle T \rangle = 608 \text{ keV/u}$







Beam acceleration test (2)

■ Entire compact-injector system




Beam acceleration test (2)

- Current of $^{12}\text{C}^{4+}$: **380 e μ A**
 - Transmission: **~96%**
- Emittance: **1.1 $\pi \cdot \text{mm} \cdot \text{mrad}$**
- Energy
 - $\langle T \rangle =$ **4.00 MeV/u**
 - $|\Delta T/T| =$ **0.3%**



Comparison of the size

HIMAC INJECTOR




Alvarez DTL

RFQ

10GHz ECR

PIG


COMPACT INJECTOR



IH-DTL RFQ


ECR

Comparison of the size




0m 5m 10m

IH-DT



10GHz ECR



Summary & conclusion

- Compact injector
 - Total length: 6 m
- Excellent performance
 - Intensity of $^{12}\text{C}^{4+}$: 380 e μA
 - Emittance : 1.1 $\pi \cdot \text{mm} \cdot \text{mrad}$
 - $|\Delta T/T| = 0.3 \%$
- First compact therapy-complex
 - Gunma University