

Recent Results from the Optical Replica Synthesizer Experiment at FLASH

Volker Ziemann, Uppsala University

for the ORS collaboration:

G. Angelova, VZ, *Uppsala University*

P. van der Meulen, P. Salén, M. Larsson, *Stockholm University*

H. Schlarb, J. Bödewadt, E. Saldin, E. Schneidmiller,

M. Yurkov, F. Lühl, A. Winter, DESY



S. Khan, *DELTA, TU Dortmund*

A. Meseck, *BESSY*


The Idea behind the ORS




- Problem: measure ultra-short bunches in the 10s of fs range: EOS, TEO, LOLA, ORS
 - too fast for electronics (10 Gs/s, 100 ps)
 - but laser folks know (autocorrelation, FROG)
- Solution: make an optical copy of the electron bunch and analyze that with laser methods.



The Seed Laser





- Er-fiber ring-oscillator (~ 1550 nm) phase locked to RF (micro-timing)
- Booster amplifier
- 2nd harmonic generation to 772 nm
- CPA 2001 regenerative amplifier on loan from Stockholm
- Pockels cell fire to let the light pulse out (macro-timing)
- 0.7 mJ/pulse, 150 fs to 2 ps
- Safety shutters (ND and other)
- Diagnostics: Frog, virtual waist






UPPSALA
UNIVERSITET



Laser Transfer Line and OS0



The Undulators




- Electromagnets
- Designed and built by Scanditronix, Vislanda
- Period 20 cm
- 5+2 periods
- 4 power supplies per magnet
- Modulator=(V)eronica
- Radiator=(H)ilda





UPPSALA
UNIVERSITET

ORS Optical Station 2









GRENOUILLE

- *Cylindrical lens* makes horizontal strip
- *Fresnel biprism* creates crossing wavefronts in thick *SHG crystal* → auto-correlator
- Effective thickness of SHG crystal varies with viewing angle → Spectrally resolved
- Second double cylindrical lens images onto camera
- Horizontally → time
- Vertically → spectrum
- GRENOUILLE USB 8-50 controlled by VideoFROG software

- Picture from Trebino's book




Experiment Preparation: Transverse Overlap




Electron orbit < 0.1 mm

Electron and Laser-
position on OTR
equal before and
after modulator




Temporal Overlap of sub-ps Electron bunch und Laser pulse

Rough adjustment on photo diode on OS1 per synchrotron radiation and laser ~ 100 ps




Fine-tuning on OS2 by observing coherent OTR of modulated electrons

- * 2nd harmonic (BG39)
- * Intensity in R.O.I



Adjust laser-timing in ~ 100 fs steps



OTR on OS2-camera while 200 fs laserpulse passes through electron bunch

PHYSICAL REVIEW SPECIAL TOPICS - ACCELERATORS AND BEAMS 11, 070702 (2008)

Observation of two-dimensional longitudinal-transverse correlations in an electron beam by laser-electron interactions

G. Angelova and V. Ziemann

Department of Physics and Astronomy, Uppsala University, 75121 Uppsala, Sweden

A. Meseck

BESSY, Albert-Einstein-Strasse 15, 12489 Berlin, Germany

P. Salén, P. van der Meulen, M. Hamberg, and M. Larsson

Department of Physics, Stockholm University, AlbaNova University Center, 10691 Stockholm, Sweden



J. Bödewadt, S. Khan, and A. Winter

Universität Hamburg, Institute of Experimental Physics, Luruper Chaussee 149, 22 761 Hamburg, Germany

H. Schlarb, F. Löh, E. Saldin, E. Schneidmiller, and M. Yurkov


DESY, Notkestraße 85, 22607 Hamburg, Germany

(Received 17 December 2007; published 24 July 2008)




090527, DIPAC


V. Ziemann: Optical Replica Synthesizer




...finally: Single-shot FROGs (preliminary, analysis in progress)



- From radiator (HILDA)
- Significant tuning with OS2 Setup
- long/short Grenouille
- First shortE/shortL because of intensity
- Here shortE/longL during SASE conditions at 700 MeV (13 nm)
- Unfortunately no simultaneous LOLA measurement
- Parasitic operation



Spin-off: ORFIR



- ORS modulator (Laser + Veronica + Chicane) generates density modulation at 780 nm
- SASE-undulator generates VUV pulse (13 nm)
- FIR-undulator emits 780 nm pulse from the same electron bunch (synchronized)
- Tested successfully in January and March in collaboration with M. Gensch's group




Conclusions

- Installed and commissioned the optical replica synthesizer in FLASH since fall 2007
- We managed to hit the electron bunch with laser
 - can be used to measure longitudinal-transverse correlations in long (few ps) bunches
- Eventually recorded online FROG traces from the short-pulse GRENOUILLE
 - unfortunately no simultaneous LOLA measurements
- Need time for parameter scans to learn system
- Electron-laser interaction spin-off: ORFIR
- ORS is scheduled to be moved to sFLASH

Backup slides







Scheme of the ORS synchronization & trigger system




Optical Station 1

Essential for timing: Laser + Synchrotron radiation




- Coupling between laser and electrons

$$\Delta U = e \int (\vec{E} \cdot d\vec{s}) = e \int E_x v_x dt$$




- Some gain, some loose, depending on initial phase





Radiator Undulator



- Electrons have longitudinal density modulation and can radiate coherently.
- Each electron slice oscillates in undulator (like an antenna) and all contributions are added in phase.
- Number of periods N determines the length of the light pulse that an electron emits → short undulator
- Need to propagate replica pulse to diagnostic section



Resolution

- 8-50: 17 fs
- 5 Periods
 - 4 $\mu\text{m}/13$ fs
- Wavefront tilt
- Dispersion in optics on OS2
- Bent(?) mirrors
- Plasma-oscillations(?)

- Grenouille: Datenblatt von Swamp Optics



GRENOUILLE model:	8-9USB	8-20USB	8-50USB	8-300USB	8-500USB
Wavelength range:	700 – 1100 nm				700 – 900 nm
Pulse-length range @ 800 nm:	~10 – ~100 fs	~20 – ~200 fs	~50 – ~500 fs	~0.3 – ~2 ps	~0.5 – ~5 ps
Pulse-length range @ 1050 nm:	~8 – ~80 fs	~15 – ~80 fs	~30 – ~100 fs	~0.1 – ~1 ps	na
Temporal resolution @ 800 nm:	3.7 fs	12 fs	17 fs	50 fs	90 fs
Temporal resolution @1050 nm:	2 fs	9 fs	13 fs	41 fs	na
Delay increment ¹ :	0.95 fs/pixel	0.85 fs/pixel	1.145 fs/pixel	11.5 fs/pixel	11.5 fs/pixel
Temporal range ³ :	336 fs	480 fs	1.9 ps	19 ps	19 ps
Spectral resolution @ 800 nm:	5 nm	4 nm	2 nm	0.23 nm	0.05 nm
Spectral resolution @1050 nm:	6.5 nm	15 nm	7 nm	0.8 nm	na
Spectral range @ 800 nm ³ :	300 nm	160 nm	50 nm	8 nm	10 nm
Spectral range @ 1050 nm ³ :	400 nm	400 nm	125 nm	20 nm	na
Pulse complexity:	Time-bandwidth product < ~10				
Intensity accuracy:	2%				
Phase accuracy:	0.01 rad (intensity-weighted phase error)				
Single-shot possible?	Call us. ²	Yes; both free-running mode & triggered single-shot are now standard.			
Sensitivity (single-shot):	Call us. ²	1 μJ			
Sensitivity (at 10 ³ pps):	500 μW (500 nJ)	100 μW (100 nJ)			
Sensitivity (at 10 ⁶ pps):	50 mW (500 pJ)	10 mW (100 pJ)			
Spatial profile accuracy:	< 0.2 % (Camera has true 8 bits and 480 x 640 pixels)				
Spatial chirp accuracy (dx/dλ):	1 μm/nm				
Pulse-front tilt accuracy (dt/dx):	0.05 fs/mm				
Required input polarization:	Any (just rotate GRENOUILLE!)				
Required input-beam diameter:	2 – 4 mm (collimated)				
Input-beam lateral-displacement tolerance:	> 1 mm				
Number of alignment knobs:	Zero				
Time to set up:	~ 10 minutes				
Dimensions (L x W x H) w/camera:	33 cm x 7.5 cm x 16.5 cm	33 cm x 7.5 cm x 16.5 cm	33 cm x 4.5 cm x 11.5 cm	45 cm x 7.5 cm x 16.5 cm	61 cm x 7.5 cm x 16.5 cm
Weight:	3 kg	3 kg	1.2 kg	3 kg	6 kg

1. At full camera resolution.

2. The Model 8-9 can be modified to allow single-shot measurement, but at a reduction in sensitivity.

3. Temporal and spectral "ranges" are the full-scale ranges, not the pulse FWHM (which is typically a factor of 2 to 3 smaller).

Comparison with LOLA (long several ps bunches)



- Simultaneous (almost, 30 min) measurement of bunch profile with transversely deflecting cavity LOLA (blue) and ORS (black).
- Initially the time calibration of LOLA was off by 20 %, now fixed.
- OD2 Neutral density filter before the Basler camera to prevent saturation
- smoothing and sqrt(ORS)
- Very good agreement of the recorded bunch length