



Science and
Technology
Facilities Council

ASTeC

Development of 6D electron beam diagnostics for novel acceleration experiments at FEBE on CLARA

Thomas Pacey,

S Mathisen, T Overton, J Wolfenden, C Swain, D Walsh, J Henderson, N Joshi, C
Tollervey, A Bainbridge, J Jones, A Pollard, E Snedden, Y Saveliev, D Angal-Kalinin

ASTeC, STFC, Daresbury Laboratory



The University of Manchester



UNIVERSITY OF
LIVERPOOL






Science and
Technology
Facilities Council

ASTeC

CLARA & FEBE

Overview

CLARA & FEBE Overview



PHASE 1:
50 MeV,
250 pC, 10 Hz

- User area Beam Area 1 (BA1)
- Beam exploitation, inc. novel acceleration

PHASE 2:
250 MeV ASSEMBLED &
INSTALLED OFFLINE

- 2022/23: Installation
- 2023/24: Commissioning

FEBE:
SHIELDED USER FACILITY

- Branch line to shielded hutch
- Flexible user exploitation space
- 100 TW class laser

PHASE 3:
100 nm FEL
NOT YET FUNDED

- Tied to decision on next stages of UK XFEL

FEBE Hutch

LASER In
Periscope


FEC #1
Laser focus
Primary IP

LASER
Out &
Dump

FEC#2
Can support flexible
array of diagnostics

Up to
600 MeV

CLARA
250 MeV



Beam parameters – “Day 1”

Parameter	High charge	Low charge
Energy [MeV]	250	250
Charge [pC]	250	5
RMS t [fs]	100	50
σ_E/E [%]	<5	<1
RMS x [μm]	100	20
RMS y [μm]	100	20
ϵ_N x [μm]	5	2
ϵ_N y [μm]	5	2

Important: all parameters to be confirmed through measurement using appropriate diagnostics.

Beam parameters – Future potential

Parameter	High charge	Low charge
Energy [MeV]	250	250
Charge [pC]	250	5 → <5
RMS t [fs]	100 → 50	50 → <50
σ_E/E [%]	<5 → 1	<1
RMS x [μm]	100 → 50	20 → 1
RMS y [μm]	100 → 50	20 → 1
ϵ_N x [μm]	5 → <5	2 → 1
ϵ_N y [μm]	5 → <1	2 → <1

R&D: Meeting these parameters will require upgraded and/or new diagnostics

Additional challenges for novel acceleration

'novel' -> plasma, laser, dielectric, THz -> high gradients, high transverse forces

- Single shot diagnostics required :
 - Beams with instability & providing feedback systems
- Meeting exacting IP parameters requires:
 - Precise measurements
 - Monitoring of shot-shot variation & longer term drifts
- Variety of experiments -> broad range of beam parameters = fundamentally challenging
- **High impact experiments need good diagnostics!**



Science and
Technology
Facilities Council


ASTeC

6D Diagnostics R&D



6D Diagnostics R&D across FEBE beam line

- Diagnostics undergoing active R&D
- All supplemented by well developed diagnostics, e.g. BPMs, YAGs, ICTs





Science and
Technology
Facilities Council



ASTeC

Longitudinal Diagnostics





Broadband PM spectrometer

- Can be deployed in second chamber for commissioning beams with energy > 600 MeV
- Energy range 50 MeV - 2 GeV
- Upstream quad. doublet + correctors for matching into dipole
- Modular design, 200mm blocks, 5 modules for 2 GeV
- Full range requires ~ 1 m scintillator screen



Dielectric wakefield streaker

- Passive, robust, single shot bunch length measurement
- RMS estimate or full profile reconstruction
- Broadband, $\sim 50 \text{ fs} > 1 \text{ ps}$
- Resolution improved using pairs orthogonally orientated plates
 - Pub. in preparation





- Resolution is given by:

$$r(s) = \frac{\sigma_Y(s)}{R_{12} \frac{dK_0(s)}{ds}}$$

$\sigma_Y(s)$ ← Transverse beam size at s
 $R_{12} \frac{dK_0(s)}{ds}$ ← Longitudinal variation in streak


Electro-optic temporal bunch profile reconstruction

- EO spectral interferometry (EOSI) system developed
 - Improved resolution over spectral decoding system (EOSD)
- Proof of concept experiments complete in CLARA BA1
 - Working across 2 – 150 pC, down to ~300fs RMS length
 - 100 pC reconstruction with crystal 1cm from beam
 - Pub. in preparation
- Future – robust fibre coupled design for FEBE experimental chamber



Coherent radiation Bunch Compression Monitors

- CTR based BCM used for BA1 for machine tuning
- Upgrade for FEBE: compact “spectrometer”
- Custom pyroelectric sensors with variable gain & active noise cancellation
- Can be replicated on arc dipoles for CSR detection
- Optional CDR target non invasive, ML driven virtual diagnostics (e.g. LPS)



Laser synchronization & TOA jitter measurement

FEBE laser oscillator will be optically synchronized to the CLARA optical clock:

- Ultralow noise optical clock (**complete**)
- Stabilized (<10 fs over 24 hr (**design and testing complete**))
- Two-colour laser-laser locking (**testing underway, all-fibre upgrade design in progress**)
- **On track to provide <10fs locking of FEBE laser to optical reference**

Beam arrival monitor (BAM) R&D:

- Developing new design for BAM with removable pick-ups
- Targeting 10fs arrival time measurement @20 pC,
- Push 10 fs @<5pC beam charge with upgraded pickups.



Science and
Technology
Facilities Council

ASTeC


Transverse Diagnostics




OTR based emittance

- Aiming for single shot emittance measurements
- Developing techniques for imaging with micro-lens array and masking with digital micromirror (DMD)
- Building up test optical systems & improving simulation tools in zemax
- Resolution limits TBC
- Collaborating with UoL and SPARC lab @ INFN

MLA method




Beamlet image




T Pacey - IBIC 2022, Krakow Poland

DMD method



DMD with mask



Angular Image

Spatial Image


R. Fiorito et al., AIP Conf. Proc. 648, (2002)
F.G. Bisesto et al, NIMA 909 (2018)

Machine Learning & Virtual Diagnostics

- Demonstrated in simulation ML prediction of IP spatial profiles as function of machine settings
- Operates in “forward” and “backward” mode
 - Arc \rightarrow IP or Dump \rightarrow IP
- Could allow for shot-shot IP prediction when IP diagnostics not available
- Future: extend to LPS prediction utilising CDR/CSR information


Forward prediction

Example optics #1



Backward prediction

Example optics #2




High dynamic range charge diagnostics


- Low charge important for some novel experiments
 - Delivers high peak current whilst maintaining transverse quality

High dynamic range

100 pC (Lowest sensitivity)




2 pC (Highest sensitivity)




Single shot <0.5pC measurement

Raw ADC waveform



Single shot after SVD



See poster contribution by S. Mathisen, Monday Session MOP32

Conclusions

- Diagnostic challenges for novel acceleration experiments
- Users will benefit from a sum of diagnostics greater than the parts
- FEBE provides strong platform to test systems for future accelerators
- We are open to collaborative opportunities for R&D

- We welcome experiments at CLARA on instrumentation
- More information on FEBE proposal call provided by 2024
- Contact directly for more information



Science and
Technology
Facilities Council

ASTeC



T Pacey - IBIC 2022, Krakow Poland

thomas.pacey@stfc.ac.uk



Science and
Technology
Facilities Council

ASTeC

Thank you