



New Developments and Status of XAIRA, the New Microfocus MX Beamline at the ALBA Synchrotron


Nahikari González

TUOAM04

- INTRODUCTION
- BEAMLINE LAYOUT
- WHITE/PINK BEAM OPTICAL COMPONENTS
 - HORIZONTAL PREFOCUSING MIRROR (HPM)
 - MONOCHROMATOR
- END STATION
 - LAYOUT OVERVIEW
 - BEAM CONDITIONING ELEMENTS
 - DIFFRACTOMETER
 - ON AXIS VIEWING SYSTEM
 - HELIUM CHAMBER
- CONCLUSIONS and NEXT STEPS

- BL06-XAIRA is the first hard X-ray microfocus beamline at the ALBA Synchrotron.
- Macromolecular crystallography (MX) beamline designed to deliver high quality data from micron-sized and/or challenging crystal samples from oscillation and fixed-target serial MX experiments.
- Beamline characteristics:
 - Beam size at sample position: $3 \times 1 \mu\text{m}^2$ (FWHM).
 - Energy Range: 3 – 15 keV (nominal 4 – 14 keV)
 - High flux ($>10^{13}$ ph/s/250mA) at 1 Å wavelength (12.4keV)
- Currently under comissioning, expecting first users in 2024.






- Windowless Optics
- Simple and stable optomechanical elements



- Maximum Flux at sample (at low energy range).
- Maximum Beam Stability at sample position.

HPM – Horizontal Prefocusing Mirror

Mirror Parameters



- Rectangular Mirror: 40x30x670mm (WxHxL)
- Horizontally deflecting mirror bender
- Bending Forces: $F_u = 438\text{N}$; $F_d = 632\text{N}$
- Incident Angle: 4,5 mrad
- Footprint (full slits aperture $0.112 \times 0.150 \text{ mrad}^2$): 2,8 x 500 mm²
- 2 adaptive optics correctors to improve slope errors

Thermal Loads

NOMINAL CASE (250mA)
Slits closed to 0.099x0.150 mrad




Total Absorbed Power

153.99 W

Peak Power Density


0.12 W/mm²

Working cond 250mA - Mirror Hrefl 4.5mrad
Total Absorbed Power = 154.0 W (Max 0.12 W/mm² at 20.100 m)




COOLING DESIGN

- Water cooling with InGa, following the procedure developed at SLAC.
- Disregarded options:
 - Internally Cooled Mirror – or Side Clamped Mirror: Mechanical interaction between cooling and mirror.
 - Trough cooling: Uncertainty in the deformation model, non-symmetric, gallinstan issues.






- Stainless Steel Tubes brazed to the long copper pads.
- Copper pads are nickel plated, to avoid copper-InGa issues
- Silicon Pads clamped to the copper pads, with Indium foils
- InGa layer between the mirror and the pads
- Thermal contact validated at the optics lab with a precision infrared camera (Optris PI640).



RESULTS



Metrology at the Optics lab




Bender Error	0,023urad rms
Radius repeatability $\Delta R/R$	0,0076%
Radius Stability (14 h)	0,0054% rms
Radius resolution (400 halfsteps)	≤ 0.0371 %
Slope Error	0.258 μ rad rms
Height Error	1.56 nm rms

Bender performance was not affected by the InGa Cooling System

Comissioning at the Beamline



Temperature increase with the front-end vertical aperture, as calculated by FEA.






H beam size measurements (edge scans at HSS pos) at different H apertures
Theoretical value achieved at narrow aperture, optimization with correctors to be done.

DESIGN CONCEPT

- Channel Cut Monochromator (CCM, Si(111)) and Double Multilayer Monochromator (DMM, Mo/B₄C) in a single mount.
- Bragg axis is 2,3mm underneath the first optical surface: the center of the beam travels along the crystals surface depending only on Bragg angle, θ .
- Due to the relatively small grazing incidence angle of the multilayers compared to the channel cut, the beam positions for the two diffracting surfaces do not overlap in the ranges:
 - 3-15 keV (θ : 41,2° - 7,2°) for the CC
 - 6-14 keV (θ : 2,27° - 0,97°) for a ML with a d spacing of 26 Å.
- Dimensions optimized to change from CC to ML just changing the Bragg rotation angle
- Minimized beam excursion (4,5mm gap – max. excursion: 2,2mm – 60 μ m at sample)



MINIMUM MECHANICS (NO TRANSLATIONS) SO AS TO MAXIMIZE BEAM STABILITY AT SAMPLE POSITION






MONO COOLING


- In-depth internal geometry optimization (CFD and FEA simulations) to make the LN₂ flow as uniform as possible and minimize turbulences while maximizing the cooling capacity.
- Multilayer cooling with LN₂ was a critical point.
 - Mo/B₄C was found to be the most suitable coating material for cryogenic temperatures
 - Cooling optimized to minimize thermal cycling on the ML substrate
- Clamping pressure was potentially an issue:
 - Essential in order to have a good thermal conductance at CC
 - It induced surface deformations on the ML optical surface



Step by step clamping procedure
measuring surface deformations with
Fizeau Interferometer





Model 6: local velocity at the center point of each micro channel



Almost constant velocity (rms) in each microchannel (~0,4m/s)

Concept validated by the first commissioning tests with beam



- Substrates changed in ~1min just rotating the Bragg axis (limited by the rotation speed)
- Horizontal beam position is almost constant ($\pm 0,1\text{mm}$).
- Measured flux $1.1 \cdot 10^{14}$ ph/s at 7,29 keV after mono (ML).



* Vertical Beam Excursion 0,34mm



Layout

Sample to detector distance **70 - 250mm**



Beam Conditioning Elements


- LOCATION: In between the KB System and the ES Helium Chamber [250mm]
- UHV Compatible, sealing with a CF16 Diamond Window ($\varnothing 2\text{mm}$ 10 μm thickness).
- Independent actuation by SmarAct PiezoGuides of the following elements:
 1. Fast Shutter [CEDRAT]
 2. 4-Blade Slits
 3. XBPM [CIVIDEC]
 4. Beam Diagnostics: Diode and Fluorescence Screen (YAG:Ce)
- XBPM vertical position is monitored by an interferometer [Qtools]
- The setup was pre-aligned at the optics lab using a telecentric lens



REQUIREMENTS


- To maximize sample position stability and minimize SoC ($\leq 0,1 \mu\text{m}$).
- To allow the positioning of the goniometer axis with respect to the beam.
- To allow the centering of the sample to rotation axis.
- To be compatible both with air and helium environments.
- To allow generating trajectories required for the experiments (i.e. helical scans, raster scans)

	Range [mm]	Resolution	Speed [mm/s]
Longitudinal (Y)	- 4 to 0,1 mm	100nm	0,5
Transversal (X)	$\pm 5\text{mm}$	100nm	0,5
Sample Centering (x&y)	$\pm 2\text{mm}$	10nm	2,5
Sample Centering (z)	$\pm 5\text{mm}$	10nm	2,5
Goniometer Rotation (Ω)	360°	0,05 mdeg	60 rpm




GONIOMETER

- Direct Drive Torque Motor; Slotless (zero cogging) from Aerotech.
- 2 Encoder Heads (Analog Encoders).
- Precision Slip-Ring (Moog) with a vibration damping coupling




1st mode 384 Hz



Temp. at goniometer 293K

SAMPLE CENTERING SETUP


- SmarActs in parallel or 3 for the vertical motion.
- Frame made in Titanium
 - Low thermal expansion coefficient.
 - Not heavy.
 - Rigid.
- Total weight: 0,6kg Easily removable. It can be exchanged by a XZ fast scanning setup (for SSX experiments).




The whole system has been calculated and optimized by FEA and CFD.

RUNOUT MEASUREMENTS

- Run out measured with a high resolution microscope.
 - At different speeds (360deg/s ; 55 deg/s; 5,5 deg/s).
 - Both with compressed air and helium.
 - With and without the slip ring.
- Motor poles can be observed.
- Behavior is similar in air and helium: concept validated.
- Run out of full setup still being optimized.




GONIO WITH AIR
50nm




GONIO WITH HE
70nm

At 360deg/s



GONIO WITH SLIP RING
65nm




FULL SETUP
242nm


NOT YET
OPTIMIZED

- Sample visualization is one of the key components of the BL.
- Based on two separate optics
- Fast exchange from one to another (< 4s)


	Description	Resolution
High Magnification Branch	High Resolution Objective	0,7 μm (fixed)
Low/Medium Magnification Branch	Parallax-free comercial system [B-ZOOM ARINAX]	1,2 μm
	\varnothing 1mm drilled hole objective	> 1,2 μm
	Splitter for 2 lens branches	



WHITE LIGHT ILLUMINATION



BLUE LIGHT ILLUMINATION




UNDER COMISSIONING

\varnothing 0.5 μm polystyrene microspheres
Scaled image showing whole FOV

- Optimized performance compatible with standard MX BL operation.
- Novel solutions for the BL optics and End Station
 - HPM: improved stability and focus
 - Dual monochromator concept: Channel Cut and Multilayer on a single mount.
 - Helium/Air compatible sample environment.
- Integrated custom design of the whole End Station.

NEXT STEPS:

- Fine commissioning of the optics: mirror correctors, interferometry system.
- Complete the installation and commissioning of the end station.



Judith Juanhuix

Damià Garriga

Isidro Crespo

Josep Nicolás

Marta Llonch

Jon Ladrera

Marcos Quispe

Carles Colldelram

Alejandro Crisol

Igors Sics

Bernat Molás

José Gabriel Centeno

Juan Luis Frieiro

Xavier Serra

Alberto Rubio

Liudmila Nikitina

Alejandro Enrique

Jose Ferrer



THANK YOU

谢谢

Q&A?

Nahikari González

TUOAM04