

Thermoacoustic Range Verification for Ion Therapy

SK Patch¹, M Kireeff-Covo², A Jackson², YM Qadadha¹,
KS Campbell², RA Albright², P Bloemhard², AP Donoghue², CR Siero²,
TL Gimpel², SM Small², BF Ninemire², MB Johnson², and L Phair²

¹UW-Milwaukee

²Lawrence Berkeley National Lab

Thermoacoustic Range Verification

inverse source problem

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar)

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ N/m}^2$$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ J/m}^3$$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ J/kg} \quad \text{kg/m}^3$$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ Gy } \rho$$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ Gy } \rho$$

$$\delta p = \Gamma \rho D \text{ dose in Gy}$$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ Gy } \rho$$

$$\delta p = \Gamma \rho D \text{ dose in Gy}$$

$\Gamma \sim 0.1$ in soft tissue

$\rho = 1000 \text{ kg/m}^3$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ Gy } \rho$$

$$\text{soft tissue } \delta p \sim 100 \text{ D } \text{dose in Gy}$$

$$\Gamma \sim 0.1 \text{ in soft tissue}$$

$$\rho = 1000 \text{ kg/m}^3$$

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ Gy } \rho$$

$$\text{soft tissue } \delta p \sim 100 \text{ D } \text{dose in Gy} \quad \Gamma \sim 0.1 \text{ in soft tissue}$$
$$\rho = 1000 \text{ kg/m}^3$$

Clinical margin: 1 mm + 3.5% target depth

Thermoacoustic Range Verification

inverse source problem

- rapid heating (stress confinement)
- pressure jump (ultrasound)
- time of flight (sonar) $v_s = 1.5 \text{ mm}/\mu\text{s}$

$$1 \text{ Pa} = 1 \text{ Gy } \rho$$

$$\text{soft tissue } \delta p \sim 100 \text{ D } \text{dose in Gy} \quad \Gamma \sim 0.1 \text{ in soft tissue}$$
$$\rho = 1000 \text{ kg/m}^3$$



Clinical margin: 1 mm + 3.5% target depth 8 mm at 20 cm

Experimental Setup

Accelerator:	88" cyclotron at LBNL 50 MeV protons, lost 1 MeV in ion chamber 2 μ A for 2 μ s on target 4.2 mm FWHM beam width 2 Gy/pulse in water target (SRIM)
Targets:	water in 112-qt LDPE container gelatin phantom mimics muscle (ultrasound)
Ultrasound Hardware:	clinical array with 96 channels, 300 μ m pitch 1000 pulses averaged

Simulations



Positive (compressional) followed by weak negative (rarefactional)



Simulations

Positive (compressional) followed by weak negative (rarefactional)



4 pC single-turn



Simulations

Positive (compressional) followed by weak negative (rarefactional)




4 pC over 2 μ s @ 10 MHz



Simulations

Positive (compressional) followed by weak negative (rarefactional)




4 pC over 2 μs @ 10 MHz



Simulations

Positive (compressional) followed by weak negative (rarefactional)




4 pC over 2 μs @ 10 MHz



Simulations

Positive (compressional) followed by weak negative (rarefactional)




4 pC over 2 μs @ 10 MHz



Simulations

Positive (compressional) followed by weak negative (rarefactional)


4 pC over 2 μ s @ 10 MHz





Simulations

Positive (compressional) followed by weak negative (rarefactional)

4 pC over 2 μs @ 10 MHz



*Obey stress confinement:
Build up pressure
faster than it runs away*






Simulations

Positive (compressional) followed by weak negative (rarefactional)


49 MeV: 4 pC in 2 μ s

4 pC in 6 μ s



Simulations


Positive (compressional) followed by
negative (rarefactional)




Simulations

Positive (compressional) followed by negative (rarefactional)

150 MeV





49 MeV



Simulations

Positive (compressional) followed by negative (rarefactional)




Henrotin, Cyclotrons 2016
Jones, Med Phys 2016

Waterbath Results

Expect 21.1 mm
range

- within 1 mm
- consistently high
- 5 realizations



1 cm left

21.7 ± 0.2

centered

21.9 ± 0.0

1 cm right

22.0 ± 0.02

Patch, et al PMB 2016


Phantom Experiment

Gelatin phantom designed to


- mimic ultrasound properties of muscle
- detect error due to gas pocket

oil range should
exceed gel range
by $\sim 10\%$

$\rho_s = 0.911$ olive oil




Phantom Results




Phantom Results

Bragg Pk in yellow; entry point in red.
Same transducer →
automatic co-registration
of image & range estimate.




oil-filled cavity




Phantom Results

Bragg Pk in yellow; entry point in red.
Same transducer →
automatic co-registration
of image & range estimate.




oil-filled cavity




Phantom Results

Bragg Pk in yellow; entry point in red.
Same transducer →
automatic co-registration
of image & range estimate.



oil-filled cavity


empty cavity



Phantom Results


Bragg Pk in yellow; entry point in red.
Same transducer →
automatic co-registration
of image & range estimate.

cavity status	distance distal (mm)	range $\mu \pm \sigma$ (mm)	N
A w/oil	64.0	18.7 ± 0.1	27
B empty	29.3	16.9 ± 0.4	28



oil-filled cavity


empty cavity



Phantom Results


Bragg Pk in yellow; entry point in red.
Same transducer →
automatic co-registration
of image & range estimate.

cavity status	distance distal (mm)	range $\mu \pm \sigma$ (mm)	N
A w/oil	64.0	18.7 ± 0.1	27
B empty	29.3	16.9 ± 0.4	28



oil-filled cavity

empty cavity




Phantom Results

Bragg Pk in yellow; entry point in red.
Same transducer →

automatic co-registration
of image & range


cavity status	distance (mm)	range $\mu \pm \sigma$ (mm)	N
A w/oil	64.0	18.7 ± 0.1	27
B empty	29.3	16.9 ± 0.4	28

10x Kellnberger, SciRep 2016
~10% difference



oil-filled cavity

empty cavity




Phantom Results

Bragg Pk in yellow; entry point in red.
Same transducer →

automatic co-registration
of image & range


cavity status	distance (mm)	range $\mu \pm \sigma$ (mm)	N
A w/oil	64.0	18.7 ± 0.1	27
B empty	29.3	16.9 ± 0.4	28

10x Kellnberger, SciRep 2016
~10% difference



oil-filled cavity

empty cavity



TA Range Ver. take-aways

- fast dose deposition increases TA signal strength and bandwidth
- off-the-shelf cardiac ultrasound array detected TA emissions from 2 Gy peak, but required overdosing the phantom
- custom ultrasound hardware sensitive to TA emissions required to detect within therapeutic dose limits

Acknowledgements

M Zolotorev & A Sessler, who coordinated first attempt @ LBNL in 2013
Dr. George Noid (MCW RadOnc) provided CT images of phantom
U.S. Department of Energy under Contract No. DE-AC02-05CH11231
Verasonics hardware: UWM Instrumentation Grant